

Table of Contents

Introduction	2
Section One	3
Installing Your Balance	3
SETUP	3
Section Two	4
Understanding Your Balance	4
KEYBOARD FUNCTIONS	4
LED INDICATORS	4
DISPLAY MESSAGES	4
Section Three	5
Using Your Balance	5
BASIC WEIGHING	5
WEIGHING WITH A CONTAINER	5
DUAL RANGE FEATURE (EL 410D & EL 4100D)	5
INTERFACING WITH A COMPUTER	6
PRINTING OUT INFORMATION	6
Appendix I	7
User Setups	7
CALIBRATION	8
SETTING FUNCTION OF MODE KEY TO PRINT OR HOLD	9
SETTING THE BAUD RATE	11
FIRMWARE VERSION	12
RESTORING THE FACTORY DEFAULT SETUPS	13
ENVIRONMENTAL CUSTOMIZATION	14
TEMPERATURE COMPENSATION CALIBRATION	16
Appendix II	19
RS-232 Serial Data Communications	19
Receiving Data Using the Immediate Print Symbol	19
The RS-232 Interface Hardware	20
Appendix III	21
Accessories and Options	21
Appendix IV	22
In Case of Difficulty	22
Appendix V	23
Specifications	23
Appendix VI	24
Limited 3 Year Warranty	24

Introduction

Thank you for purchasing a Setra precision balance. The fine workmanship and durable construction should provide years of reliable service. While your balance is easy to operate, it is advisable to read this guide carefully before use. It is designed to help you perform weighing and related operations quickly and accurately.

This manual is divided into four major sections. Section One, "Installing Your Balance," explains where to put your balance, how to level it and install the pan, and how to get started. Section Two, "Understanding Your Balance," explains the various keys, displays, and messages you will encounter while using your balance. Section Three, "Using Your Balance," provides the detailed instructions necessary to perform various operations. Following Section Three are appendices which include initial balance setups, RS-232 information, accessory information, troubleshooting, specifications and a warranty statement.

Typographical conventions used in this manual include the following:

1. **BOLD UPPER CASE CHARACTERS** indicate specific keys on the balance keyboard.
2. "Quotation marks" enclose messages seen on the balance display.

Installing Your Balance

SETUP

1. Locating your balance...

You should treat your balance as you would any piece of precision equipment, locating it on a clean, dry surface; away from extremes of draft and vibration.

2. Leveling your balance...

The bubble level is located underneath the loading pan. Adjust the two front feet until the bubble is centered in the circle.

3. Installing the weighing pan...

Next, place the pan support over the shaft in the center of the unit. Then place the pan over the pan support.

4. Connecting your balance to an AC outlet...

To power up your balance, insert the appropriate end of the AC adaptor into an electrical outlet. Now insert the connector end into the back of the balance. The balance draws very little current from an outlet, and should always remain plugged in. The display may be turned off by holding down the **TARE** key for five seconds. This keeps your balance always ready to use with no "warm up" time.

5. Display test...

Each time your balance is powered on, it automatically runs through a display test lasting about thirty seconds, showing all possible segments which might be displayed, and then setting itself to zero. The balance is then ready to be used. For best accuracy allow the balance to warm up for 5 minutes prior to using or calibrating.

Understanding Your Balance

KEYBOARD FUNCTIONS

TARE/POWER Assigns the pan and whatever is currently being weighed a value of zero. To turn the display off, hold this key down for five seconds. To turn the display on press this key again.

MODE Factory programmed to function as a print key, this key can be changed to access different print modes or to activate the hold feature.

LED INDICATORS

ok Reading shown is stable.

g Reading shown is given in grams.

DISPLAY MESSAGES

----- The balance is developing a stable reading.

UnAbLE The balance is unable to perform your requested operation. Press the **TARE** key and select another operation.

HHHHHH The weight on the pan exceeds the capacity of the balance.

LLLLLL The pan is not properly seated or has been removed.

Using Your Balance

BASIC WEIGHING

To weigh a sample on your balance, use the following procedure:

1. Press the **TARE** key to zero the display.
2. Place the object(s) to be weighed on the pan.
3. Wait for the “ok” indicator, then read the weight from the display.

WEIGHING WITH A CONTAINER

To weigh objects or liquids without including the weight of the container, use the following procedure:

1. Place the empty container on the pan. Press the **TARE** key, the scale will display “-----” and return to zero.
2. Wait for the “ok” indicator. Place or pour objects or liquids into the container.
3. Wait for the “ok” indicator, the net weight will be displayed.

DUAL RANGE FEATURE (EL 410D & EL 4100D)

Each time the **TARE** key is pressed the display will present the higher of the two possible weight resolutions (0.001g for the low capacity balance, and 0.01g for the high capacity balance). The dual range feature provides this higher resolution up to 100 grams for the low capacity model, and up to 1000 grams for the high capacity model. Whenever these weights are exceeded, the balance will automatically reduce the resolution by one order of magnitude. The higher resolution can be recovered by either removing some weight, or by pressing the **TARE** key and resetting the display to zero. Regardless of the weight on the pan the higher resolution can always be displayed by pressing **TARE**. Pressing the **TARE** key does not increase the capacity of the balance, thus the total weight on the pan must always be equal to or less than the rated capacity.

INTERFACING WITH A COMPUTER

Your balance has a male DB9 RS-232 serial port and is designed to interface with computer equipment. If your balance is connected to a computer, follow the instructions in Appendix II.

PRINTING OUT INFORMATION


Your balance is designed to print out the displayed weight when connected to an optional serial printer. To print using the CoStar SETRA 250 thermal receipt/label printer, follow the instructions below:

1. Connect the printer's AC adaptor to the proper electrical outlet.
2. Make sure the printer is turned on (as indicated by the printer's green light). If the printer is not on, press the blue power button on the front of the printer.
3. Load the appropriate paper or label stock into the printer. (See printer instructions for details).
4. Connect the printer to the balance's RS-232 connector using the cable provided.
5. Perform the necessary weighing procedures on the balance.
6. Press the **MODE** key on the balance.

NOTE: The **MODE** key must be set up to perform the print function. This is the factory default. If using label stock, the form feed command must be programmed in User Setups (see Appendix I). When using a printer other than the CoStar SETRA 250, set the baud rate and parity of your balance to match the printer (see Appendix I, User Setups to select the print mode, baud rate and parity).

User Setups

The outline below represents the balance's menu structure. To enter the menu press and hold the **TARE** key and then press the **MODE** key. To view the current menu options, press the **TARE** key repeatedly. To select the displayed option, press the **MODE** key.



CALIBRATION


To perform a span calibration, use the following procedure:

STEP 1.



Press the **TARE** key to zero the balance.

STEP 2.



Hold down the **TARE** key then press the **MODE** key. Release both keys, the display will read “SPAN” for span calibration.

STEP 3.



Press the **MODE** key to enter the span calibration menu. The balance’s full capacity will be displayed. You can calibrate the balance at full capacity or at a lower calibration point. To view the set calibration points, press the **TARE** key until the desired weight is displayed.

To restore the factory calibration, press the **TARE** key to display “UnSPAN” then press the **MODE** key. To escape, press the **TARE** key to display “ESCAPE” then press the **MODE** key.

STEP 4.




Place on the pan calibration weight(s) equal to the displayed weight and press the **MODE** key. The balance will display “ACAL” and then return to the normal weighing mode. If the display reads “nOCAL”, check the weights and weighing environment. This calibration procedure can only correct for $\pm 1\%$ span shift.

SETTING FUNCTION OF MODE KEY TO PRINT OR HOLD

The **MODE** key may be set to either print or hold a displayed weight. The selectable print functions are: “stable print” which will only print once a stable reading is attained, “interval print” which may be programmed to print at predetermined time intervals, “instant print” which will print immediately after the **MODE** key is pressed (note: the reading may not be stable). The number of line feeds may also be set for label printing. The print function is set separately from the line feed setup, i.e., set the print function first then re-enter the print menu to program the number of line feeds. The **MODE** key can also be set to the hold feature in which the balance attains a stable weight and locks it in the display until a key is pressed.

To set the function of the **MODE** key, use the following procedures:

STEP 1:



Hold down the **TARE** key then press the **MODE** key. Release both keys, the display will read “SPAN”.

STEP 2.



Press the **TARE** key again and the display will read “FUNC” for function.

NOTE: To escape anytime during this procedure, press the **TARE** key until “ESCAPE” is displayed, then press the **MODE** key.

STEP 3.



Press the **MODE** key to enter the function menu. The display will read “PrInT”. Proceed to step 4a. to program the **MODE** key to perform a print command. Go to step 4b. to program the **MODE** key to perform the hold function.

STEP 4:



a. Press the **MODE** key to enter the print menu. The display will read "STABLE".

For stable print

Perform step 4a. then press the **MODE** key to select the stable print mode. The balance will return to the normal weighing mode.

For instant print

Perform step 4a. then press the **TARE** key once to display "InSTAn" for instant and then press the **MODE** key. The balance will return to the normal weighing mode.

For interval print

Perform step 4a.; press the **TARE** key twice to display "InTer" for interval print and then press the **MODE** key. Proceed to step 5.

For line feed

Perform step 4a. then press the **TARE** key three times to display "LnEFD" and then press the **MODE** key. Proceed to step 6.



b. To select the hold feature press the **TARE** key to display "HOLD" then press the **MODE** key as shown in step 4a. The balance will return to the normal weighing mode. Each time the **MODE** key is pressed the weight will lock on the display until a key is pressed.

STEP 5:



To view the predetermined print intervals (in seconds) press the **TARE** key repeatedly. When the desired time interval is displayed, press the **MODE** key. (Select zero for continuous printing.) The balance will then return to the normal weighing mode. Pressing the **MODE** key will print the displayed weight after each selected time interval (e.g., every 90 seconds). To interrupt the interval printing press the **MODE** key again. To re-activate, press the **MODE** key.

Note: Print intervals can vary up to $\pm .2$ seconds depending on weight variations.

STEP 6:




To view the preset number of line feeds available (1-18) press the **TARE** key repeatedly. To program a form feed command for the CoStar SETRA 250 printer, select "LF 0" as the number of line feeds. When the desired number of line feeds is displayed, press the **MODE** key. The balance will then return to the normal weighing mode.

SETTING THE BAUD RATE

The balance is capable of interfacing with a wide variety of computer devices. To set the baud rate (the rate at which the scale communicates with a computer/printer) and parity, use the following procedure:

STEP 1.



Hold down the **TARE** key then press the **MODE** key. Release both keys, to display "SPAN".

STEP 2.



Press the **TARE** key twice, the display will read "bAUD" for baud rate.

Note: To escape anytime during this procedure, press the **TARE** key until "ESCAPE" is displayed and press the **MODE** key.

STEP 3.



Press the **MODE** key to enter the baud rate menu. The display will read 300. To view the other baud rates press the **TARE** key repeatedly.

STEP 4.



When the desired baud rate is displayed, press the **MODE** key to select it. The display will then read "PARITY".

STEP 5.



Press the **MODE** key to enter the parity menu. The display will read "nONE" for no parity. To view the parity menu press the **TARE** key.

STEP 6.




When the desired parity (none, odd, even) is displayed, press the **MODE** key. The balance will then return to the normal weighing mode.

FIRMWARE VERSION

The operating software in your balance has a reference number. To display this number, follow the procedure outlined below.

STEP 1.



Hold down the **TARE** key then press the **MODE** key. Release both keys, the display will read "SPAn".

STEP 2.



Press the **TARE** key repeatedly to display "rEF" for reference number.

Note: To escape anytime during this procedure, press the **TARE** key until "ESCAPE" is displayed and press the **MODE** key.

STEP 3.



Press the **MODE** key to display the software reference number.

STEP 4.




Press the **TARE** key to return to the normal weighing mode.

RESTORING THE FACTORY DEFAULT SETUPS

The many features described in this section allow the user to customize the balance to suit a particular application. However, in doing this it is also possible to inadvertently set up the balance in such a way that it does not operate as expected. To reset the factory defaults to: **MODE** key will print a stable reading and then advance the printer to the next line, auto response rate, 2400 baud, no parity, follow the steps below.

Note: Restoring the factory defaults will return your balance to all of the factory span and temperature calibration settings. You **MUST** recalibrate (span) your balance after restoring the factory defaults. If you are experiencing a temperature induced offset, you should also run the temperature compensation calibration procedure.

STEP 1.



Hold down the **TARE** key then press the **MODE** key. Release both keys, the display will read "SPAN".

STEP 2.



Press the **TARE** key repeatedly until the display reads "FACdEF" for factory defaults.

Note: To escape anytime during this procedure, press the **TARE** key until "ESCAPE" is displayed and press the **MODE** key.

STEP 3.



Press the **MODE** key to restore the set ups to the original factory defaults. The balance will display "buSY" and then return to the normal weighing mode.

ENVIRONMENTAL CUSTOMIZATION

Most high precision balances are used in stable environments, away from drafts and vibration. You may customize the balance with regards to its environmental conditions. The balance's response rate can be set to fast, slow or auto, depending on your application.


The balance is set to the auto mode from the factory. This mode is recommended for most applications and will provide good stability and response.

The fast mode will provide a faster response than the auto mode and is suitable for applications such as filling to a target weight. In environments where draft and vibration are present, the displayed weight may not be as stable as required.

The slow mode will provide maximum stability especially in a noisy or drafty environment; however, the response rate may be slower.

To set the display response rate for your application, perform the following steps:

STEP 1.



Hold down the **TARE** key and then press the **MODE** key. Release both keys, the display will read “SPAN”.

STEP 2.



Press the **TARE** key repeatedly until the display reads “CUSTON” for the customization menu.

Note: To escape anytime during this procedure, press the **TARE** key until “ESCAPE” is displayed and press the **MODE** key.

STEP 3.



Press the **MODE** key to enter the customization menu and the display will read “SPEEd” for the balance’s response rate setting.

STEP 4.



Press the **MODE** key to enter the response rate menu and the display will read “FAST” for a fast response rate.

For fast response

Perform Step 4., then press the **MODE** key to select the fast response rate. The balance will return to the normal weighing mode.

For slow response

Perform Step 4., then press the **TARE** key once to display “SLO”, then press the **MODE** key to select the slow response rate. The balance will return to the normal weighing mode.

For auto response

Perform Step 4., then press the **TARE** key twice to display “AUTO”, then press the **MODE** key to select the auto response rate. The balance will return to the normal weighing mode.

TEMPERATURE COMPENSATION CALIBRATION

If your balance is frequently moved from one location to another or if the room temperature varies more than a few degrees during the day, the balance may display a slight temperature induced offset. Your Setra balance is capable of compensating for this offset so that varying temperatures will not affect its performance.

If your balance is used in different locations, leave it where the temperature is either higher or lower than the operating environment for at least four hours. Otherwise, perform this calibration first thing in the morning or last thing in the evening to expose the balance to the largest possible temperature fluctuation. Then, place the balance in its operating environment and run this calibration procedure. Over the next three hours it will monitor any weight fluctuations that occur and record these temperature induced variations into memory. During this calibration period the balance cannot be used and must not be disturbed. Once this calibration procedure is completed, the balance will be able to internally compensate for the offset.

To perform the temperature compensation calibration:

STEP 1.


Leave the balance for at least four hours in the initial environment as explained above. The balance does not have to be powered on.

STEP 2.

Move the balance to its operating environment. Unplug the balance (if powered on) for 10 seconds and then power on. The balance will count down and should display 0.00.

Note: The 0.00 reading may not be very stable if the balance is warming up or cooling down.

STEP 3.



Remove the pan and pan support. Hold down the **TARE** key then press the **MODE** key. Release both keys, the display will read "SPAn".

STEP 4.



Press the **TARE** key repeatedly until the display reads "CUSTON" for the customization menu.

Note: To escape anytime during this procedure, press the **TARE** key until "ESCAPE" is displayed and then press the **MODE** key.

STEP 5.



Press the **MODE** key to enter the customization menu. The display will read "SPEEd" for the balance's response rate setting.

STEP 6.



Press the **TARE** key once and the display will read "CAL T" for the temperature calibration menu.

STEP 7.



Press the **MODE** key to select the temperature calibration menu and the display will read "FAC TP" for the factory temperature calibration. Proceed

to Step 8 to calibrate the balance to its operating environment. Proceed to Step 10 to restore the factory temperature calibration.

STEP 8.



Press the **TARE** key to display “USER TP”, then press the **MODE** key. The display will read “HOTCAL” to begin temperature calibration. Let the balance sit, without being disturbed, for approximately three hours to monitor the temperature variations and weight fluctuations.

STEP 9.



At the end of approximately three hours, the balance will display a number such as “0.03”. This number represents the maximum temperature induced offset that occurred as the balance acclimated itself to the operating environment. Press the **TARE** key and the balance will return to the normal weighing mode.

STEP 10.



To restore to the factory temperature calibration, press the **MODE** key. The balance will restore the factory temperature calibration and return to the normal weighing mode.

RS-232 Serial Data Communications

The balance keyboard functions can be accessed via the RS-232 interface. The following commands are available:

t = tare function
p = print function
= immediate print (see below)

Receiving Data Using the Immediate Print Symbol

When a balance is connected to a computer, it is suggested that immediate print (“#”) be used. In response to this command the balance will transmit whatever number or message appears on the balance display. The “string format” output is shown below:

+/- 1 2 3 4 5 6 . c0 c1 c2 c3 CR LF

The first six digits represent the number field. A sign (+ or -) always precedes the number and a decimal point is always transmitted. Numbers less than six digits long are preceded by spaces. (Messages, when transmitted, are sent in the number field.)

NOTE: The position of the decimal point will depend on the readability and units the balance is displaying. The sign will be adjacent to the leading digit.

“c0” is a space.

“c1” is also a blank space as long as the balance is in the automatic display response mode. If the display response rate has been changed to fast, the “c1” will be an “F”, and if the display response rate has been changed to slow, then “c1” will be an “S”.

“c2” is the “units” character. It describes the units of the number being transmitted. Your balance will transmit “G” for grams.

“c3” is the “stability” character. This character corresponds to the “ok” indicator on the display. A (space) means the reading is not stable. “S” means the reading is stable.

The immediate print output is always transmitted with a carriage return and line feed. If the balance is set to a specific number of line feeds, these will be transmitted with a carriage return.

The RS-232 Interface Hardware

Although Setra balances can communicate with almost any RS-232 device, the built-in interface does not include the complete protocol. Only the transmit and receive lines of the standard interface are used. This should not present any interfacing problems in most applications.

The data format is:

- 1 start bit
- 8 data bits including parity
- 1 stop bit
- 10 bits per frame (framing errors ignored)

Note: The balance will transmit using the parity selected; however it does not check the parity it receives. Use an RS-232 cable to connect the external device to the balance, or construct one following the instructions below.

Connect a high quality, shielded cable with a DB9S (D-Subminiature 9 pin female connector) using the following pinout:

1 2 3 4 5	PIN	DESCRIPTION
• • • • •	2	TXD - scale transmits data
• • • •	3	RXD - scale receives data
6 7 8 9	5	GRD - signal ground

Note: "Handshake" signals, such as "Clear To Send" (CTS), are not used. The peripheral must have a minimum buffer (15 characters).

Shielded cables must be used with this unit to ensure compliance with the Class A FCC limits.

Computers which require handshaking need a connection between two pins on the computer's connector named DTR and DSR (Data Terminal Ready and Data Set Ready). CTS may also need to be jumpered to RTS at your computer interface (Clear To Send and Request To Send). The maximum recommended cable length is 15 meters. The cable can be longer if it has < 2,500 pF capacitance. The load impedance of the device connected should be between 3,000 and 7,000 ohms with no more than 2,500 pF shunt capacitance.

Accessories and Options

The following accessories are available for the EL series of precision balances. Contact your dealer for current price information.

Part #	Description
401160	CoStar SETRA 250 Printer - 2.25" wide receipt/label thermal printer. Direct connect for weight only printing. Dimensions: 4"W x 8"L x 6.25"H (100mm x 200mm x 160mm).
401907	Receipt Paper - 300' Roll (100 m).
401908	Labels - 2.25"w x 1.25"H (57mm x 32mm). Adhesive backed labels, 1000/roll.
407910	3 Way Sliding Glass Door Draftshield - Detachable glass and stainless steel compartment provides static-free protection from draft-created instability. Three sliding doors, including a top opening for tall vessels makes it easy to fill to a desired weight. For use on 200 and 410 gram balances only.
407210	Replacement Glass Cylinder Draftshield – Detachable glass draftshield provides static-free protection from draft-created instability. Removable stainless steel cover (p/n 407212-02) has a small opening for taller vessels. For use on 200 and 410 gram balances only.
407900	Security Lock Kit - Secure your balance to a bench or table top. Combination lock (not provided) can be used to permit removal of instruments.
407901	Dustcover - Chemically resistant rubber membrane cover protects balance from powder dust and spills.

In Case of Difficulty

If the balance will no longer follow your instructions, unplug it from its power source; then, plug the unit in again. If any unusual messages appear during warm up, or if the balance does not return to normal operation, contact your dealer or Setra Systems, Inc.

If the balance displays “-----” for an extended period of time, or the displayed reading is unstable, too much vibration or draft may be present. Relocate the balance away from the source of vibration or shield the balance from draft. If it continues, service may be required.

If the balance displays “nOCAL” during calibration, check to make sure you are using the correct calibration weight. (This calibration procedure can only correct for $\pm 1\%$ span shift.)

If you are experiencing difficulty in calibrating or printing, check the set up of the balance. To restore the factory default setups:

1. Hold down the **TARE** key then press the **MODE** key. Release both keys, the display will read "SPAn".
2. Press the **TARE** key repeatedly until the display reads "FACdEF" for factory defaults.
3. Press the **MODE** key. The balance should display “buSY” and return to the normal weighing mode.

If the RS-232 interface does not function correctly, first make certain the RS-232 cable is securely fastened to both the balance and the peripheral device. Next, reset the balance's baud rate and parity to match the external device and make certain the data formats are identical. If data transmission or reception is still not possible, check that the cable is the correct type. (It may be necessary to “cross” the receive and transmit lines of the interface. That is, the receive line of the balance must connect to the transmit line of the external device. Likewise, the transmit line of the balance must connect to the receive line of the external device. Special cables can be purchased for this purpose from a computer dealer.) See Appendix II for detailed information.

If the problem still persists, contact your dealer or Setra Systems, Inc.

Specifications

Model	EL-200S	EL-2000S	EL-410S	EL-4100S
Capacity (g)	200	2000	410	4100
Readability (g)	.001	.01	.001	.01
Repeatability (g)	± .001	± .01	± .001	± .01
Linearity (g)	± .002	± .02	± .002	± .02
Keyboard	two key	two key	two key	two key
Pan Size	4" dia.	6.25" dia.	4" dia.	6.25" dia.
Draftshield	standard	N/A	standard	N/A
Security Lock Kit	option	option	option	option
Size WxDxH (in.)	7 x 5.75 x 11	7 x 3 x 11	7 x 5.75 x 11	7 x 3 x 11
Weighing Units	grams	grams	grams	grams
Interface	Bidirectional RS-232, (standard on all models)			
Display	0.57" LED (all models)			
Power	6-9 VDC			

Model	EL-410D	EL-4100D
Capacity (g)	100/410	1000/4100
Readability (g)	.001/.01	.01/.1
Repeatability (g)	± .001/.01	± .01/.1
Linearity (g)	± .002/.01	± .02/.1
Keyboard	two key	two key
Pan Size	4" dia.	6.25" dia.
Draftshield	standard	N/A
Security Lock Kit	option	option
Size WxDxH (in.)	7 x 5.75 x 11	7 x 3 x 11
Weighing Units	grams	grams
Interface	Bidirectional RS-232, (standard on all models)	
Display	0.57" LED (all models)	
Power	6-9 VDC	

Span Range:	± 1% of factory calibration
Display Update Speed:	Up to 5 times per second
Selectable Baud Rates:	300 600 1200 2400 4800 9600
Selectable Parity:	none/odd/even
Pan:	Stainless Steel
Operating Temperature:	40° to 100° F (5° to 43° C)
Shipping Weight:	10 lbs.

Limited 3 Year Warranty

Setra Systems, Inc. warrants the EL balances it manufactures to be free from defects in material and workmanship. Upon return, transportation charges prepaid, to an Authorized Service Center within three (3) years of the date of purchase, Setra or its Authorized Agent will repair or replace, at its option, any balance which it determines to contain defective material or workmanship and will return said balance to purchaser, transportation prepaid. Setra shall not be obligated, however, to repair or replace balances which have been repaired by unauthorized parties, abused, improperly installed, altered, or otherwise misused or damaged, even if by accident, in any way. Setra will not be responsible for any dismantling, reassembly or reinstallation charges.

Nothing in this warranty shall be construed as a warranty for merchantability or fitness for any specific use or purpose, and this warranty is in lieu of all other warranties, express or implied. Setra shall not be held liable under the terms of this warranty for any special, indirect, incidental or consequential damages claimed in connection with the balances' performance or availability.

Setra Systems, Inc.
Weighing Systems Division
159 Swanson Road • Boxborough, MA 01719

FCC WARNING

Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area may cause harmful interference in which case the user will be required to correct the interference at their own expense.

Shielded cables must be used with this unit to ensure compliance with the Class A FCC limits.

CE Compliance Statement

Top Loading Balances

Setra Systems, Inc. manufactures high accuracy weighing equipment and distributes its products worldwide. In an effort to meet the demands of the market, Setra Systems continues to improve its products following the guidelines set forth by various governmental bodies and regulatory agencies. One such guideline has been established by the European Community, effective 1 January 1996. This guideline is known as the EMC Directive.

The EMC Directive requires compliance with various published standards for electromagnetic compatibility. Such standards specify testing to be done in the following areas:

1. EN50081-1 and EN50081-2 1994
Electromagnetic Compatibility, Generic Emission Standard
Industrial, Scientific and Medical Equipment
 - a) EN55011, Radiated Radio Frequency (RF)

2. EN50082-1 1995
Electromagnetic Compatibility, Generic Immunity Standard
Residential, Commercial, Light Industrial Environments
 - a) IEC801-2, Electro Static Discharge (ESD)
 - b) IEC801-3, Radiated Radio Frequency (RF)
 - c) IEC801-4, Transient Burst

Setra Systems, Inc. has tested and certified its top loading precision balances for compliance with the EMC Directive. This certification is evidenced by the CE label on the product and the Declaration of Conformity on file. There may be applications where EMC effects are not a concern. However, some applications may require optimum performance under certain conditions specified in the standards. The results of testing verifies that this product complies with the EMC guidelines, notwithstanding the fact that there may be some effect on accuracy under certain conditions. This Application Note contains the results of the EMC Directive testing for this product to verify compliance with applicable standards.

Any questions on this Application Note can be directed to the supplier of this product in your country. See the following page for the results of the EMC testing.

Results of Compliance Testing for EMC Directives

Products: Top Loading Balances

Model EL-200S
 Model EL-410S
 Model EL-2000S
 Model EL-4100S
 Model EL-410D
 Model EL-4100D

EMC Directive 1994 EN 50081-1 and EN 50081-2 Electromagnetic Compatibility
 Generic Emission Standard for Industrial, Scientific and Medical Equipment

Test No.	Level	Criteria	Results
EN55011 Radiated RF	30-230 MHz	< 30 dB uV/m	*Pass
EN55011 Radiated RF	230-1000 MHz	< 37 dB uV/m	*Pass

EMC Directive EN 50082-1 1995 Electromagnetic Compatibility
 Generic Immunity Standard for Residential, Commercial, Light Industry

Test No.: IEC801-2 ESD	Level: 8kV Air Discharge
Criteria: B - Must Resume Normal Operation After Test	Results: Pass
Test No.: IEC801-2 ESD	Level: 4kV Direct Discharge
Criteria: B - Must Resume Normal Operation After Test	Results: Pass
Test No.: IEC801-3 Radiated RF	Level: 3V/m 27-500 MHz
Criteria: A - Must Operate As Specified	Results: **Pass
Test No.: IEC801-4 Transient Burst	Level: 5 ms Burst, 5KHz
Criteria: B - Must Resume Normal Operation After Test	Results: Pass

* Use of ferrite bead on power lead required.

** Typical susceptibility to radiated frequencies is < 0.0005% FS. At certain discreet frequencies there is the potential for a maximum effect of 0.18% FS.